

## Product Datasheet | HLA-C\*07:02 easYmer® kit

<b>Application</b>	<p>easYmers® are peptide loadable HLA class I monomers.</p> <p>easYmers® can be used to generate peptide-HLA monomers with your choice of peptide. The monomers can easily be tetramerized with fluorophore conjugated streptavidin and used to analyse T cells by flow cytometry. The easYmer reagent can also be used to evaluate specific peptide-HLA I binding.</p>			
<b>Catalog no.</b>	1122-01			
<b>Size</b>	<b>Test</b>	<b>Vol.</b>	<b>Folding Buffer</b>	<b>Control Peptide</b>
<input type="checkbox"/>	20	40 µl	1 vial (0.5 ml)	1 vial
<input type="checkbox"/>	50	100 µl	1 vial (1 ml)	1 vial
<input type="checkbox"/>	150	300 µl	1 vial (1 ml)	2 vials
<input type="checkbox"/>	500	1000 µl	1 vial (3 ml)	3 vials
<b>Allotype</b>	<p>HLA-C*07:02</p> <p>peptide receptive, biotinylated in Tris/Maleate pH 7, 30% Glycerol</p>			
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.			
<b>Shelf life</b>	<p>easYmers® are stable at -20°C for 12 months.</p> <p>See expiration on vial</p>			
<b>Folding buffer</b>	Tris/Maleate pH 7			
<b>Control Peptide</b>	<p>FRAPNTREL (20 nmol lyophilized peptide per vial )</p> <p>HLA-C*07:02 binder and positive control for evaluation of peptide-HLA folding. Dissolve each peptide vial in 20 µl DMSO resulting in a peptide conc. of 1 mM. Store at -20°C</p>			
<b>Peptide source</b>				

**For Research Use Only (RUO)**