

Product Datasheet | HLA-B*40:01 easYmer® kit

Application	easYmers® are peptide loadable HLA class I monomers.				
	easYmers® can be used to generate peptide-HLA monomers with your choice of peptide. The monomers can easily be tetramerized with fluorophore conjugated streptavidin and used to analyse T cells by flow cytometry. The easYmer reagent can also be used to evaluate specific peptide-HLA I binding.				
Catalog no.	1082-01				
Size	Test	Vol.	Folding Buffer	Control Peptide	
	<input type="checkbox"/> 20	40 µl	1 vial (0.5 ml)	1 vial	
	<input type="checkbox"/> 50	100 µl	1 vial (1 ml)	1 vial	
	<input type="checkbox"/> 150	300 µl	1 vial (1 ml)	2 vials	
	<input type="checkbox"/> 500	1000 µl	1 vial (3 ml)	3 vials	
Allotype	HLA-B*40:01 peptide receptive, biotinylated in Tris/Maleate pH 7, 30% Glycerol				
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.				
Shelf life	easYmers® are stable at -20°C for 12 months. See expiration on vial				
Folding buffer	Tris/Maleate pH 7				
Control Peptide	GEAMDTISV (20 nmol lyophilized peptide per vial)				
	HLA-B*40:01 binder and positive control for evaluation of peptide-HLA folding. Dissolve each peptide vial in 20 µl DMSO resulting in a peptide conc.of 1 mM. Store at -20°C				
Peptide source	YFV NS4A 21-29				

For Research Use Only (RUO)