

Product Datasheet | HLA-A*11:01 easYmer® kit

Application	<p>easYmers® are peptide loadable HLA class I monomers.</p> <p>easYmers® can be used to generate peptide-HLA monomers with your choice of peptide. The monomers can easily be tetramerized with fluorophore conjugated streptavidin and used to analyse T cells by flow cytometry. The easYmer reagent can also be used to evaluate specific peptide-HLA I binding.</p>			
Catalog no.	1018-01			
Size	Test	Vol.	Folding Buffer	Control Peptide
	<input type="checkbox"/> 20	40 µl	1 vial (0.5 ml)	1 vial
	<input type="checkbox"/> 50	100 µl	1 vial (1 ml)	1 vial
	<input type="checkbox"/> 150	300 µl	1 vial (1 ml)	2 vials
	<input type="checkbox"/> 500	1000 µl	1 vial (3 ml)	3 vials
Allotype	<p>HLA-A*11:01</p> <p>peptide receptive, biotinylated in Tris/Maleate pH 7, 30% Glycerol</p>			
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.			
Shelf life	<p>easYmers® are stable at -20°C for 12 months.</p> <p>See expiration on vial</p>			
Folding buffer	Tris/Maleate pH 7			
Control Peptide	<p>ALNTITNLK (20 nmol lyophilized peptide per vial)</p> <p>HLA-A*11:01 binder and positive control for evaluation of peptide-HLA folding. Dissolve each peptide vial in 20 µl DMSO resulting in a peptide conc. of 1 mM. Store at -20°C</p>			
Peptide source				

For Research Use Only (RUO)